

IABP Handover Tool

- IABP has at least 2 points of securement. Must be sutured proximally. Statlock permissible as distal securement.
- IABP Size Matches Adapter Tubing
- Operating Mode, Trigger, and Balloon Inflation (augmentation) match provider orders
- Console plugged into red outlet and battery light indicates full charge
- Helium tubing free of blood, flecks, or condensation
- Central lumen pressure bag is at least three feet above patient and fully inflated
- Central lumen is transduced with 2 units/ml heparinized saline
- Central lumen transducer is leveled and zeroed to phlebostatic axis
- Central lumen tubing has occlusive redcaps on all ports
- Extra length is removed from transducer tubing (maximum 8ft of tubing)

Please refer to the *IABP SOP* for more information regarding IABP Standards of Care pertaining to assessments, frequency, interventions, and management:



Fiberoptic Catheter?
Scan Here:



Last Updated: 11/9/2022