

Open Heart Surgery Room Set Up

- 1-EKG cables
- 1-Red dots package
- 1-SpO2 cable with box
- 1-BP cable
- 1-A-line cable
- 1-CVP cable with box
- 1-PA cable with box
- 1-CO cable with box
- 1-1L D5NS
- 2-500mL NS
- 1-Primary (non-pump) IV tubing
- 2-3 way stopcock (only need 1 if packaged together)
- 1-10mL syringe (for cardiac output set up)
- 1- Red top lab tube (for cardiac output temperature sensor-make sure tub has fluid in it)
- 1-Vamp
- 1-Transducer holder
- 1-Foley stat lock
- 1-60mL syringe (for OG)
- 1-Primary pump tubing (for medline)
- 2-Suction caniters (If ETT is a hi-lo need 3)
- 2-Suction regulators (if ETT is a hi-lo need 3)
- 1-Mouth care kits
- 1-Long suction tubing (for CT)
- 1-Graduated cylinder (for emptying foley bag)

*Restraints are in blue cart if patient needs them.

*If patient is cold will need a bear hugger with blanket.